


Experiment 1

Analysis of Plasmid DNA by Restriction Digestion and Agarose Gel Electrophoresis


This experiment will provide an introduction to restriction enzyme digestion and gel electrophoresis, the two most fundamental techniques in recombinant DNA methodology. You will learn how to cut DNA with a restriction endonuclease and analyze the products on a gel. At the end of the experiment you should be able to construct a restriction map of a plasmid.

The plasmid pUC9 (pages 41 and 42) is 2.7 kb long and contains the gene for β -lactamase (*bla*), which confers resistance to ampicillin, the *lac* repressor, operator, and promoter, and the alpha peptide (containing the amino-terminal end) of the β -galactosidase (*lacZ*) gene. When in the appropriate host strain, one that contains the carboxyl end of *lacZ*, pUC9 produces functional β -galactosidase. Within the *lacZ* gene is a closely-spaced series of restriction sites, called the multiple cloning site, which may be used as insertion points for foreign DNA fragments. Insertion into *lacZ* disrupts the β -galactosidase and cells carrying this plasmid, therefore are lac^- .

When a DNA fragment is inserted into a plasmid, it can be found in either of two orientations. In the diagram below, for example, the narrow line represents a plasmid into which a DNA fragment has been inserted into restriction site A. The plasmid and the insert each have a site for enzyme B. On the insert, B is asymmetric. In the left orientation, the two B sites are close to one another (proximal) and in the right orientation, they are far apart (distal).



If both of these plasmids are cut with enzyme A, the vector and insert fragments will be regenerated. Both orientations will produce exactly the same fragments. However, if the plasmids are cut with enzyme B, the resulting fragments will be different. The proximal orientation will produce a large and a small fragment while the distal orientation will produce two moderately sized fragments:



pRHR30 is a derivative of pUC9 in which a fragment of human DNA has been inserted into the *EcoRI* site in *lacZ*. The fragment has an asymmetric *HindIII* site which is either proximal or distal to the *HindIII* site on pUC9. In this experiment, you will cut pUC9 and pRHR30 with *EcoRI*, *HindIII*, and *BamHI* and construct a plasmid map. The results of this experiment will be combined with Experiment 2 for your first formal lab report.

Note that in addition to the restriction digests done in this experiment, you will also perform several restriction digests that will be used for Experiment 3.

Procedure

1. Prepare the series of restriction digests described below according to the procedure on page 12. You will use 5 μ l of DNA and 15 μ l of reaction mix.



Tube	DNA	Enzyme	Buffer (2 μ l)	Water	* Dye
1	pUC9 (3 μ l)	none	none	17 μ l	Orange
2	“	<i>HindIII</i>	Med	15 μ l	“
3	“	<i>BamHI</i>	NEB <i>BamHI</i>	“	“
4	“	<i>EcoRI</i>	High	“	“
5	“	<i>PstI</i>	NEB 3	“	“
6	1 kb ladder (3 μ l)	none	none	17 μ l	“
7	λ <i>HindIII</i> (3 μ l)	none	none	17 μ l	Blue
8	pRHR30 (5 μ l)	none	none	15 μ l	“
9	“	<i>HindIII</i>	Med	13 μ l	“
10	“	<i>BamHI</i>	NEB <i>BamHI</i>	“	“
11	“	<i>EcoRI</i>	High	“	“
12	“	<i>PstI</i>	NEB 3	“	“

*Tracking dye (5 μ l) is added **AFTER** the incubation period

2. Prepare the restriction digests and heat inactivations for steps 1 and 2 of Experiment 3.
3. Place samples 4, 5, 6, 8, 9, and 10 in a 37^o water bath for 40 minutes.
4. While you are waiting for your DNA's to digest, pour a 1% agarose gel according to the procedure on page 17.
5. At the end of 40 minutes, remove the samples and add 5 μ l of tracking dye to each sample.
6. Load the entire 25 μ l onto the gel and run at 120 volts for about 90 minutes.
7. Stain your gel with ethidium bromide and photograph.

Analysis

1. Measure the distance that each band migrated from the origin and record the results in a table.
2. Graph the migration of your molecular weight standards (samples 1, 2, 11, & 12) on semi-log paper. Plot $\log kb$ versus distance. This will be your standard curve. I will examine your standard curves to make sure that they are plotted correctly.
3. When your standard curve is correct, plot the distances of your cut samples on the curve and determine the size of each fragment. Record the migration distances and molecular weights in a table.
4. Construct a map for pRHR30



LIST OF M13mp VECTORS

M13mp1	1 2 3 4 5 6 THR MET ILE THR ASP SER ATG ACC ATG ATT ACG GAT TCA
M13mp2	1 2 3 4 5 6 THR MET ILE THR ASN SER ATG ACC ATG ATT ACG AAT TCA EcoR I
M13mp7/pUC7	1 2 3 4 5 (1 2 3 4 5 6 7 8 9 10 11 12 13 14) 6 THR MET ILE THR ASN ser pro asp pro ser thr cys arg ser thr asp pro gly asn SER ATG ACC ATG ATT ACG AAT TCC CCG GAT CCG TCG ACC TGC AGG TCG ACG GAT CCG GGG AAT TCA EcoR I BamH I Sal I Pst I Sal I BamH I EcoR I Acc I Hinc II Acc I Hinc II
M13mp8/pUC8	1 2 3 4 5 6 (1 2 3 4 5 6 7 8 9 10 11) 7 8 THR MET ILE THR ASN SER arg gly ser val asp leu gln pro ser leu ala LEU ALA ATG ACC ATG ATT ACG AAT TCC CCG GGA TCC GTC GAC CTG CAG CCA AGC TTG GCA CTG GCC EcoR I BamH I Pst I Hind III Sma I Xma I Sal I Acc I Hinc II
M13mp9/pUC9	1 2 3 4 (1 2 3 4 5 6 7 8 9 10 11) 5 6 7 8 THR MET ILE THR pro ser leu ala ala gly arg arg ile pro gly ASN SER LEU ALA ATG ACC ATG ATT ACG CCA AGC TTG GCT GCA GGT CGA CGG ATC CCC GGG AAT TCA CTG GCC Hind III Pst I Sal I BamH I Sma I EcoR I Acc I Xma I Hinc II
M13mp10/pUC12	1 2 3 4 5 6 (1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16) 7 8 THR MET ILE THR ASN SER ser ser pro gly cys arg pro leu glu ser thr cys ser pro ser leu ala LEU ALA ACC ATG ATT ACG AAT TCG AGC TCG CCC GGG GAT CCT CTA GAG TCG ACC TGC AGC CCA AGC TTG GCA CTG GCC EcoR I Sst I BamH I Xba I Sal I Pst I Hind III Sma I Xma I Acc I Hinc II
M13mp11/pUC13	1 2 3 4 (1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16) 5 6 7 8 THR MET ILE THR pro ser leu gly cys arg ser thr leu glu asp pro arg ala ser ASN SER LEU ALA ACC ATG ATT ACG CCA AGC TTG GGC TGC AGG TCG ACT CTA GAG GAT CCC CGG GCG AGC TCG AAT TCA CTG GCC Hind III Pst I Xba I BamH I Sst I EcoR I Sal I Sma I Xma I Acc I Hinc II
M13mp18/pUC18	1 2 3 4 5 6 (1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18) 7 8 THR MET ILE THR ASN SER ser ser val pro gly asp pro leu glu ser thr cys arg his ala ser leu ala LEU ALA ATG ACC ATG ATT ACG AAT TCG AGC TCG GTA CCC GGG GAT CCT CTA GAG TCG ACC TGC AGG CAT GCA AGC TTG GCA CTG GCC EcoR I Sst I Kpn I BamH I Xba I Sal I Pst I Sph I Hind III Sma I Xma I Acc I Hinc II
M13mp19/pUC19	1 2 3 4 (1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18) 5 6 7 8 THR MET ILE THR pro ser leu his ala cys arg ser thr leu glu asp pro arg val pro ser ser ASN SER LEU ALA ATG ACC ATG ATT ACG CCA AGC TTG CAT GCC TGC AGG TCG ACT CTA GAG GAT CCC CGG GTA CCG AGC TCG AAT TCA CTG GCC Hind III Sph I Pst I Xba I BamH I Kpn I Sst I EcoR I Sal I Sma I Xma I Acc I Hinc II